Technical Compliance Sheet for "Optical Tweezers Facility"

Please fill-up the compliance sheet as per the given instructions completely and accurately as partially filled compliance sheet shall not be accepted for evaluation of technical bids. Any modification in the format of the compliance sheet is not permitted and mandatory cell should not be left blank. Partially/inaccurately filled compliance sheet may lead to rejection of technical bid.

S. No.	Technical Specification	Complied (YES/NO/EQUIVALENT /BETTER)	Explanation and details of supporting document(s) It is mandatory to give details. If you have filled YES, please give the exact specifications provided by the company, and a link/reference to the company brochure. If you have filled EQUIVALENT or BETTER, please articulate the <i>how</i> with appropriate company specification details and links to their documents. You can attach the company document/s supporting the technical specifications as Annexure/s.
1	The system should consist of a dual-beam optical tweezers with option of magnetic twisters fully integrated with multimode imaging system (with multiple possible imaging modes such as confocal fluorescence/fluorescence/TIRF/Widefield/IRM/DIC etc.), optical table with vibration isolation, along with required computers, electronics, software, accessories required for operation, data collection and data analysis.		

	All technical specifications of the multimode imaging	
	system and optical table with vibration isolation being	
	offered by bidder should be clearly stated without any	
	ambiguity. Offered multimode imaging system should	
	be compatible with electronics, software, accessories	
	required for operation, data collection and data	
	analysis of the optical tweezers facility.	
1A	Fully motorized and computer-controlled inverted	
	fluorescence research microscope with multiple port	
	design:	
i	Optics for bright field, fluorescence and DIC	
	observation	
ii	It is preferred to have a motorized universal condenser	
	turret with NA 0.55 or better; motorized 6 (or better)	
	position objective turret with 4 or more DIC slots;	
	motorized filter turret with at least 6 positions for	
	band pass fluorescent filters for sample visualization.	
	Band pass filters for GFP/FITC, DAPI, TRITC/	
	Rhodamine/ cy3 should be quoted.	
iii	Binocular tube for sample visualization.	
iv	Compatible C mount adapter for multiple ports (at	
	least two for simultaneous attachment of two	
	cameras).	
v	CMOS USB 3.0 b/w camera for C-mount	
vi	The ports should be camera-ready with all required	
	optics for epifluorescence/bright-field imaging/DIC	
	imaging. Further, it should be possible to send 100%	
	light to each port	
vii	Built-in motorized z-focus drive capable of step size of	
	10 nm or better without encoder, and step size of 20	
	nm or better with encoder.	
viii	High power LED for diascopic illumination.	
	Epifluorescence adaptor with a remote-	
	controlled/computer-controlled fluorescence	

	excitation source/mercury light source. Extra lamps	
	(one for each type) should be quoted for both.	
ix	Plano apochromatic spring-loaded objectives: 10x/0.45	
	or better NA (air), 20x/0.75 or better NA (air),	
	40x/0.75or better (air), 60x/1.40NA or better (oil	
	immersion), 100x 1.45NA or better (oil immersion)	
	with complete DIC accessories for all objectives. All	
	objectives should be corrected for UV, visible and IR	
	wavelengths.	
x	High resolution plano apochromatic confocal and	
	multiphoton grade spring-loaded objectives: 60x/1.27	
	or better NA (water immersion) suitable for bright field	
	imaging, and 60x/1.27 or better NA (oil immersion) for	
	super-resolution imaging.	
xi	Additional 40X (NA 1.1 or better) water-immersion	
	confocal grade apochromatic objective dedicated for	
	FCS/FCCS should be quoted.	
xii	Additional 60X (NA 1.49 or better) oil-immersion and	
	100X (NA 1.49 or better) oil-immersion confocal grade	
	apochromatic objective dedicated for TIRF should be	
	quoted.	
1B	Motorized XYZ Stage:	
i	Motorized and programmable X-Y specimen scanning	
-	stage including controller with universal sample holder	
	for slides, Petri dish and multi-well plate.	
ii	XY motorized stage with 15x15 mm movement, High	
	precision XY stage with minimum 100x100 micron	
	range of motion with 10 nm or less resolution in X-Y	
	with required controllers.	
iii	Piezo/Galvo Z-stage (mountable on the X-Y stage) with	
	minimum 100 micron range of motion, 10 nm or less	
	resolution in Z, 5 nm repeatability, 10 nm step size,	

	and 0.5 to 12.5 mm/s speed range.	
iv	Hardware-enabled auto-Z drift correction focus drift	
	compensation device for long-term time-lapse imaging applications	
1C	Temperature and CO2-controlled incubation stage for live cell imaging:	
i	Programmable/software-controlled fully automatic on- stage incubation system with temperature, 95% humidity and 5% CO2 control (capable of operating with 100% CO2 gas supply). Large incubation chamber for temperature control and a smaller chamber for CO2 control. Active CO2 and temperature control with sensors. The chambers should be compatible with operation with the piezo/galvo stage.	
ii	The environment chamber and stage should accommodate standard chamber slides, slides with cover-glass, multiple round bottom cover-glass dishes of various sizes and multi-well plates. The system should support live cell imaging up to 72 hours.	
1D	Sample Holder for optical Tweezer applications:	
i	Sample holder for slides, Petri dish or custom fluid cells	
ii	Motorized XY sample positioning (minimum 15x15mm)	
iii	XYZ Piezo nanopositioning sample scanner. 100 micron x 100 micron x 100 micron	
iv	Multichannel laminar flow cell with temperature Control	

v	Temperature range from room temperature 42°C	
vi	Sealed dual cover slip design	
vii	Up to 5 parallel flow channels fed through tube connectors	
viii	Starter kit for laminar flow cell	
ix	Dual syringe pump for fluid flow	
x	Fully integrated software control	
1E	High-performance active vibration isolation optical table:	
i	Active 1800 mm x 1200 mm floating anti-vibration optical table with minimum thickness 305 mm for the complete confocal and optical tweezers system with at least two active and two precision tuned passive dampers	
ii	Working surface 4.8 mm thick 400 series ferromagnetic stainless steel with Surface Flatness± 0.1 mm over 600 mm square, Mounting HolesM6, Mounting Hole Pattern25 mm grid	
iii	Equipped with Pneumatic Vibration Isolator with Automatic Re-leveling feature of minimum 28 in. Height along with required compatible accessories such as air-compressor etc.	
1F	Lasers (including control electronics, shutter and manipulator for fiber output):	

i	Diode laser with 488 nm	
ii	DPSS laser with 561 nm	
iii	Diode laser with 640 nm	
iv	It is also preferred to have a Diode laser with 405 nm	
V	A maximum deviation of +/- 2 nm from above- mentioned wavelength is acceptable.	
vi	Output power at the fiber end for each laser: 15 mW or better with built-in Acousto-Optic Tunable Filter (AOTF)	
vii	All visible lasers should be connected to the scan head through fibre optic cable and should be computer- controlled through AOTF for fast laser switching and attenuation in synchronization with the laser scanner.	
1G	Laser scan head:	
i	Confocal laser point scanning unit should be quoted. It should be capable of line-scanning in 2D (space-time) as well.	
ii	Scanner should have laser ports to connect to lasers and should have either low angle of incidence dichroic for high efficient excitation laser suppression and for applications like FCS/FCCS.	
iii	Standard filter for quoted wavelengths i.e. 405, 488, 561, 640 (with +/- 2 nm deviation as mentioned above)	
iv	Motorised and computer-controlled continuously- variable confocal pinhole with software control.	
V	Scanner: galvano scanners, Pixel size: at least 4096 x 4096 pixels, 2 fps or better (512 x 512 pixels, bi- direction, 1x zoom)	

vi	Zoom: ~1-1000x continuously variable	
vii	Scan mode: X-Y, X-T, X-Z, XY rotation, Free line, Line-Z.	
viii	The laser scanner should have dual scan capability of fast volumetric regions bleaching/activation/ablation & normal scan for Imaging to conduct experiments like FRAP, FRET, FLIP, photo activation/conversion.	
ix	It should also be capable of real ROI-based scanning for fast time-lapse experiments.	
1H	Detectors:	
i	The detection unit should have dual detection capability with intensity-based confocal imaging as well as spectral confocal imaging.	
ii	Should have minimum 3 detectors (of which at least 2 should be highly sensitive GaAsP/HyD) or a 32-channel high sensitive spectral detector or equivalent with minimum 45% quantum efficiency. The system should be capable of simultaneous detection and separation of at least 3 fluorophores with built-in confocal detectors (at least 2 fluorophore imaging with high sensitive GaAsP/HyD). All the detectors should be actively cooled for best SNR image over extended period of time with excellent dynamic range for low light imaging applications.	
iii	It is preferred to have minimum 2-channel (violet/green, green/red and CFP/YFP combinations) NDD GaAsP/APD(SPAD)/HyD for multiphoton detection in reflected fluorescence light path should be quoted.	
iv	The spectral dispersion of the emission light should be of latest technology with high efficient separation. The system should be capable of online separation of	

	autofluorescence and bleed through.	
v	Capability for lambda-scan should be quoted.	
	The system should be capable of recording emission spectra with minimum spectral resolution of 5nm of better.	
vi	CCD/EMCCD camera from a leading brand compatible with confocal microscope should be quoted.	
1I	Computer with monitor, keyboard and mouse:	
	The below specifications are the minimum requirements. Windows 10 Pro 64bit, CPU-Intel Xeon W-2125 (4.0GHz, quad core, 8.25 MB, 2666 MHz) or higher, RAM- 64GB or higher, HDD-1st HP Z Turbo G2 512GB PCIe M.2 SSD or better, 2nd SATA HDD 2TB or better, Optical Drive-Super Multi drive, up to x16 speed or higher, Extension slot-Two PCI Express 3.0 (x16) slots (one slot to be used for graphics), One PCI Express 3.0 (x8) slot, Two PCI Express 2.0 (x4) slot, LAN port 10/100/1000 Network/Interface x 2 (for connection to controller, for connection to external LAN), High resolution large LCD/LED screen with 30" diagonal, 2560x1600 pixel resolution, dual monitor configuration recommended Computer should be from a reputed company factory and not of assembled type.	
Dual-	Beam Optical Tweezers:	
2	The system should be optimized for ultra-sensitive single-molecule experiments, trapping and tracking of nano-particles.	
3	Laser module: Continuous wave (no time-sharing)	

	tweezers, wavelength of 1064 nm, 3W laser power or	
	higher, lasers with different output power preferred	
4	One of two beams with adjustable power distribution	
	preferred.	
5	Force detection, independent on both optical traps	
	with stability preferred.	
6	Steering of beams using advanced piezo mirror or	
	AOD. System with AOD option for beam steering	
	preferred. Technical specification piezo mirror/AOD	
	should be clearly stated.	
i	It is preferred that two independent steering units	
	should be configured at least One with standard Piezo	
	and one with AOD to provide maximum flexibility.	
ii	De-Multiplexing of detector signals for up to 8 traps for	
	quasi simultaneous force measurements with AODs.	
iii	Decoupled two-channel detection [XY (quadrant	
	photodiode) and Z (photodiode)] for optimized Z	
	sensitivity.	
iv	Non-linearity of trap stiffness and sensitivity <5% over	
	the entire field of view.	
V	3.5 MHz bandwidth InGaAs photodiodes optimal for	
	use with 1064 nm laser source.	
vi	Motorized discrete (4 step) dynamic range adjustment	
7	The system should include a force-detection unit,	
	allowing force measurements. System with lowest	
	force resolution preferred. Minimum force resolution	
	to be clearly stated.	
8	Fully integrated software suite for system operation	
	with possibility of scripting	
9	All possible applications of offered system should be	
	clearly stated. Cost should include	
	Installation/commission and initial operator training.	
	On-site application training and application start-up	
	kits/consumables must be supplied by bidder.	

10	Cost should include 2 years on site community	
10	Cost should include 3 year on-site comprehensive	
	warranty.	
11	Cost of Annual maintenance contract (AMC) for	
	additional 2 years subsequent to the warranty period	
	should be quoted. Bidder should also include full	
	description of services to be provided during the AMC	
	period.	
12 a	Bidders should also enclose:	
	List of similar equipment installations made in India or	
	abroad during last 5 years	
12 b	Bidders should also enclose:	
	Internationally published papers using the offered	
	model should be provided. Enclose complete citations	
	of the published papers.	
Magne	etic twister:	
_		
i	Supply and integration of a Magnetic Twister to apply	
	torque and rotate magnetically sensitive particles is	
	desirable.	
ii	It should be able generate homogeneous magnetic	
	field of few hundred mT at distance of 1mm.	
iii	Magnetic twister should provide precise rotational	
	control with accuracy of at least 4 degree using	
	stepper motors and should have integrated software	
	control.	
iv	Magnetic twister should not disturb bright field and	
	fluorescence imaging.	
FCS/F	CCS Imaging System:	
i	System should be capable of performing FCS/FCCS	
	(auto and cross correlation) with high sensitive	
	detectors.	
ii	This can be quoted as either card/RICS-based or dual-	

	channel GaAsP/APD-based.	
iii	There should be minimum after-pulsing for the complete visible range in live cells and solutions for a wide range of dyes and proteins like FITC, R6G, GFP, Cy3, Cy5, R6G/Cy3, GFP/RFP, Cy3/Cy5, CFP/YFP etc.	
iv	All laser lines for confocal imaging should be capable of working in FCS/FCCS mode.	
v	Measurement software for auto and cross-correlation capabilities should be quoted.	
vi	Additional software modules for diffusion with models for rotational, flow, free as well as anomalous translational diffusion should be provided.	
vii	Software for enhanced fitting with algorithms for global and interactive fit.	
viii	Software for photon counting histogram should also be quoted.	